

TASS Management



Nick Mamalis, MD
Salt Lake City, Utah

"ASICO offers the most comprehensive report on management of TASS to date."

- Nick Mamalis, MD, University of Utah, Salt Lake City, Utah


Designs for TASS Management

AE7-0029 Cleaning I/A Handle




Works on all popular phaco consoles.


Works with various threaded I/A tips.



Skinny I/A



Ball I/A




Completely disassembles for a thorough cleaning.

Recommended tips include:

Skinny I/A: AE7-3050, AE7-3051, AE7-3052, AE7-4050, AE7-4051

Ball I/A: AE7-3060, AE7-3061, AE7-3062, AE7-3060S, AE7-3061S, AE7-3062S


AE-9063CSP Royale III Cleaning Injector



One hand plunger style for controlled release.

Ideal for the Monarch® C and D cartridges.

*Trademark of Alcon Laboratories



Completely disassembles for a thorough cleaning.

US Patent 7,047,649